



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Ratiolab centrifuge tube 15 ml, conical, sterilized

Special Price

€123.00 was

€175.38

**plus VAT &
Shipping**

Product Images



Description

Centrifuge tubes made of polypropylene, with conical bottom or with standing rim, in sterilized design. The sterile centrifuge tubes are certified RNase, DNase, endotoxin, pyrogen and ATP-free. The centrifuge tubes are provided with red screw cap, printed graduation and title block. Temperature resistance -80°C to 121°C

Additional Information

No.	RB-0044
Manufacturer (Brand)	Ratiolab
EAN	4058072984410
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	15 ml
sterile	Yes
Type closure/ lid	Screw cap
TBST MAX	121 °C
TBST MIN	-80 °C
Floor outside	conical
centrifugable up to	7600 g

