



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® PE snap ring caps ND 11 blue, 6 mm
hole, septum 1.0 mm natural rubber red-
orange/TEF, 100 pcs./pack**

**€17.20
plus VAT &
Shipping**

Product Images



Description

Blue snap ring caps ND11 made of PE, with hole 6 mm.

With septum natural rubber/TEF

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness: 1.0 mm

Additional Information

No.	7-0752
Manufacturer (Brand)	neochrom
EAN	4058072030032
Transport temperature	Room temperature
Color	blue
Type closure/ lid	Snap lock
Closures and lids features	Reference to untested autosampler compatibility in shopping cart with mounted septum suitable for autosampler suitable for multiple injections cap with septum
Quantity components	two
DM Hole	6 mm
DM Closure	11 mm
Thickness septum	1 mm
Color closure/ lid	blue
Strength cap	hard (HDPE)
Hardness	60 shore A
Material gasket / septum	Natural rubber
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

