



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® snap ring caps PE blue, 6mm hole  
ND11, Red Rubber/PTFE beige Agilent quality  
45°shore A, 1.0 mm thickness, 100 pcs./pack**

Special Price

**€13.00** was

**€19.50**

**plus VAT &  
Shipping**

## Product Images

## Description

---

Blue snap ring caps ND11 made of PE, with hole 6 mm.

With septum Red Rubber/PTFE

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness: 1.0 mm

## Additional Information

No.	7-0754
Manufacturer (Brand)	neochrom
EAN	4058072030056
Transport temperature	Room temperature
Color	blue
Type closure/ lid	Snap lock
Closures and lids features	Reference to non-tested autosampler compatibility in shopping cart with mounted septum suitable for autosampler cap with septum septum with IH quality
DM Hole	6 mm
Quantity components	two
Hardness	45 shore A
Material gasket / septum	Polytetrafluoroethylene (PTFE)
Closure / lid properties	with hole with mounted septum
Color closure/ lid	blue
DM Closure	11 mm
Thickness septum	1 mm
Strength cap	soft (LDPE)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	110 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

