



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® snap ring caps ND11, PE red with  
hole, silicone white/PTFE red, 100 pcs./pack**

**€33.90  
plus VAT &  
Shipping**

## Product Images



## Description

---

Red snap ring caps ND11 made of PE, with hole 6 mm.

With septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single use injections

Thickness: 1.0 mm, septum silicone white/PTFE red 1.3 mm

## Additional Information

No.	7-0757
Manufacturer (Brand)	neochrom
EAN	4058072030087
Transport temperature	Room temperature
Color	red
Type closure/ lid	Snap lock
Closures and lids features	Reference to non-tested autosampler compatibility in shopping cart with mounted septum with high purity septum for critical analyses suitable for autosampler cap with septum
Quantity components	two
DM Hole	6 mm
DM Closure	11 mm
Thickness septum	1.3 mm
Color closure/ lid	red
Strength cap	hard (HDPE)
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

