



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® PE snap ring caps ND11 blue,
septum silicone white/PTFE blue, cross-slit
hardness 55° shore A, 100 pcs./pack**

**€46.50
plus VAT &
Shipping**

Product Images



Description

Blue snap ring caps ND11 made of PE, with hole 6 mm.

With septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single use injections

Thickness: 1.0 mm, septum silicone white/PTFE red 1.3 mm

Additional Information

No.	7-0759
Manufacturer (Brand)	neochrom
EAN	4058072030100
Transport temperature	Room temperature
Color	blue
Type closure/ lid	Snap lock
Closures and lids features	Reference to non-tested autosampler compatibility in shopping cart with mounted septum with slotted septum for blunt. thick needles suitable for autosampler cap with septum
DM Hole	6 mm
Quantity components	two
Hardness	55 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Closure / lid properties	with hole with mounted septum
Color closure/ lid	blue
DM Closure	11 mm
Thickness septum	1 mm
Strength cap	hard (HDPE)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

