



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lid for 2-3845

€15.90
plus VAT &
Shipping

Product Images



Description

Lid made of crystal-clear polystyrene for instrument trays

Temperature-resistant up to approx. +70°C.

Order No. 2-3847

Suitable for tray 2-3845

Accessories for melamine instrument trays

Additional Information

No.	2-3847
Manufacturer (Brand)	neoLab
EAN	4058072001513
Transport temperature	Room temperature
Color	transparent (clear)
Material	Polystyrene (PS)
Wide	290 mm
Height	30 mm
Length	160 mm
Depth	160 mm
Type closure/ lid	Overlay closure
Closures and lids features	with handle
Closure / lid properties	with handle / knob
for part no.	710430

