



neoLab® Lid for 2-3855

Product Images



Description

Clear polystyrene lid for instrument trays

Temperature resistant up to approx. +70°C.

Order No. 2-3858

Suitable for tray 2-3855

Accessories for melamine instrument trays

Additional Information

No.	2-3858
Manufacturer (Brand)	neoLab
EAN	4058072001605
Color	transparent (klar)
Material	Polystyrol (PS)
Wide	210 mm
Height	20 mm
Length	110 mm
Depth	110 mm
Type closure/ lid	Aufliegender Verschluss
Closure / lid properties	geschlossen
for part no.	714082

