




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Nalgene replacement lid with quick release, for suction device large

**€244.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Quick release cap made of PP

Suitable as vacuum drop system for balloon bottles.

Replacement lid, with 2 quick release connectors (1/4") for 6 mm hose, suitable for 4 and 5 Ltr bottle.

## Additional Information

---

No.	CD-0429
Manufacturer (Brand)	Nalgene
old neoLab article no.	2-9341
EAN	4058072219086
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	White
Material	Polypropylene (PP)
sterile	No
Type closure/ lid	Screw cap with olives Screw cap with quick release connector
Closure / lid properties	with gasket
for DM hose inside 1	6.0 mm
DM Closure	120 mm
Color closure/ lid	white

