



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® silicone lid, size L, pink

**€3.77**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### DURAN® silicone lids

- Can be used for a wide range of vessels and their different shapes for covering and sealing
- The autoclavable DURAN® silicone lids close flush (precisely fitting) with the tube and thus protect e.g. against dust ingress or splashing of the reaction solution.
- The different colors and the additional text field on the flap allow samples to be visibly and clearly marked.
- Easy to clean, dishwasher safe and reusable

## Additional Information

---

No.	SD-2917
Manufacturer (Brand)	DURAN
EAN	4058072365073
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	83,7 mm
Color	pink
Material	Silicone
sterile	No
Type closure/ lid	Slip closure size L
suitable for	Bottles. round beakers
Texture	smooth
Strength cap	flexible soft (silicone)
TBST MAX	180 °C
TBST MIN	-40 °C
for bottle shape	around

