



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® silicone lid, size S, cyan

€2.29
plus VAT &
Shipping

Product Images



Description

DURAN® silicone lids

- Can be used for a wide range of vessels and their different shapes for covering and sealing
- The autoclavable DURAN® silicone lids close flush (precisely fitting) with the tube and thus protect e.g. against dust ingress or splashing of the reaction solution.
- The different colors and the additional text field on the flap allow samples to be visibly and clearly can be visibly and clearly marked.
- Easy to clean, dishwasher safe and reusable

Additional Information

No.	SD-2919
Manufacturer (Brand)	DURAN
EAN	4058072365097
Transport temperature	Room temperature
autoclavable at 121°C	Yes
diameter	42,2 mm
Color	Cyan
Material	Silicone
sterile	No
Type closure/ lid	Slip closure size S
suitable for	Bottles. round beakers
TBST MAX	180 °C
TBST MIN	-40 °C
Texture	smooth
Strength cap	flexible soft (silicone)
for bottle shape	around

