



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® plastic lid, gray-white, suitable for 2-3183

**Price on
request
plus VAT &
Shipping**

Product Images



Description

Plastic lid in signal white RAL 9003. The lid is equipped with stainless steel hinges. Easy to remove for cleaning.

Order No. 2-3185: Suitable for trash bag stand Order No. 2-3183

Additional Information

No.	2-3185
Manufacturer (Brand)	neoLab
EAN	4058072001025
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	gray
Material	Polypropylene (PP)
sterile	No
Color closure/ lid	gray
for part no.	468639

