



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Elma plastic lid for ultrasonic equipment

€37.00
plus VAT &
Shipping

Product Images



Description

Robust ultrasonic equipment for efficient cleaning of glass parts, metal parts, cutlery and plastic items.

- Stainless steel oscillating tank
- Ultrasonic power transducer 37 kHz
- optimal sound field distribution due to switchable frequency modulation
- Display of set and remaining time
- LED temperature display
- Timer and continuous operation
- FTZ-tested

Order No. 2-4146

suitable for ultrasonic unit 2-4148.

Additional Information

No.	GH-0016
Manufacturer (Brand)	Elma
old neoLab article no.	2-4146
EAN	4058072260859
Transport temperature	Room temperature
Color	blue
Wide	240 mm
Height	25 mm
Length	300 mm
Depth	300 mm
for part no.	821098

