



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® closing caps made of PE, with star,
for PP microvials, 100 pcs./pack**

**€9.90
plus VAT &
Shipping**

Product Images

Description

Closure caps for micro test tubes can be put on and removed without tools. They fit all tubes, including most standard vials.

Material: PE

With star for easier piercing.

Additional Information

No.	2-1011
Manufacturer (Brand)	neochrom
EAN	4058072072209
Color	transparent (milky)
Type closure/ lid	Push-on/plug-in closure
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart suitable for autosampler
Closure / lid properties	pierceable (with star) closed
Color closure/ lid	transparent (milky)
Material closure/ lid	Low-density polyethylene (LDPE)
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

