



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® PE snap ring caps ND11 white,
septum silicone white/PTFE blue, cross-slit
hardness 55° shore A, 100 pcs./pack**

**€46.50
plus VAT &
Shipping**

Product Images



Description

White snap ring caps ND11 made of PE, with hole 6 mm.

With septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single use injections

Thickness: 1.0 mm, septum silicone white/PTFE red 1.3 mm

Additional Information

No.	7-0758
Manufacturer (Brand)	neochrom
EAN	4058072030094
Transport temperature	Room temperature
Color	transparent (milky)
Type closure/ lid	Snap lock
Closures and lids features	Reference to non-tested autosampler compatibility in shopping cart with mounted septum with slotted septum for blunt. thick needles suitable for autosampler cap with septum
Quantity components	two
DM Hole	6 mm
DM Closure	11 mm
Thickness septum	1 mm
Color closure/ lid	transparent (milky)
Strength cap	hard (HDPE)
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

