



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® PE caps ND11 w. Hole, white, hard,  
silicone white/PTFE blue, cross-connected, 100  
pcs./pack**

**€46.50  
plus VAT &  
Shipping**

## Product Images



## Description

---

### **Colored snap ring caps ND11 made of PE, with hole 6 mm.**

1. with septum natural rubber/TEF and Red Rubber/PTFE

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness: 1.0 mm

2. with septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single injections

Thickness: 1.0 mm, septum silicone white/PTFE red 1.3 mm

## Additional Information

No.	EC-1323
Manufacturer (Brand)	neochrom
EAN	4058072086084
Transport temperature	Room temperature
Color	transparent (milky)
Closures and lids features	with mounted septum with slotted septum for blunt. thick needles suitable for autosampler closure with septum
DM Hole	6 mm
Quantity components	two
Hardness	55 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
DM Closure	11 mm
Thickness septum	1 mm
Color closure/ lid	transparent (milky)
Strength cap	soft (LDPE)
Material closure/ lid	High-density polyethylene (HDPE)
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

