



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lid without handle (PC) for neoStackRack, white

€5.36
plus VAT &
Shipping

Product Images



Description

Lid for neoStapelRack®

Without handle

Material: PC white, autoclavable

Temperature resistant from -130°C to +130°C

Additional Information

No.	2-3001
Manufacturer (Brand)	neoLabLine
EAN	4058072000134
Transport temperature	Room temperature
Color	White
Material	Polycarbonate (PC)
Wide	150 mm
Height	11 mm
Length	146 mm
Depth	146 mm
Type closure/ lid	Overlay closure
Closure / lid properties	closed
for part no.	401435

