



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Lid with handle, for neoStackRack, blue

€7.01
plus VAT &
Shipping

Product Images



Description

Lid for neoStapelRack®

With handle

Material: PC blue, autoclavable

Temperature resistant from -130°C to +130°C

Additional Information

No.	2-3003
Manufacturer (Brand)	neoLabLine
EAN	4058072000158
Transport temperature	Room temperature
Color	blue
Material	Polycarbonate (PC)
Wide	150 mm
Height	15 mm
Length	146 mm
Depth	146 mm
Type closure/ lid	Overlay closure
Closures and lids features	with handle
Closure / lid properties	with handle / knob
for part no.	401435

