



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® PE Push-On Cap ND11 with hole, natural rubber red-or.TEF transp., 100 pcs./pack


Special Price

€9.00 was

€13.50

plus VAT &
Shipping

Product Images



Description

For vials ND11 with rolled rim.

Transparent PE caps 13 x 7.5 mm with 4.5 mm hole,
Septum natural rubber red-orange/TEF transparent 1.0 mm thick.

Additional Information

No.	7-0653
Manufacturer (Brand)	neochrom
EAN	4058072029180
Transport temperature	Room temperature
Color	transparent (clear)
Type closure/ lid	Push-on/plug-in closure
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart suitable for autosampler cap with septum
DM Hole	4.5 mm
Hardness	60 shore A
Material gasket / septum	Natural rubber Polytetrafluoroethylene (PTFE)
Closure / lid properties	with mounted septum with hole
Color closure/ lid	transparent (clear)
DM Closure	11 mm
Thickness septum	1 mm
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

