



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® PE Push-On cap 8 mm natural
rubber red-orange/TEF transparent 100
pcs./pack**

**€18.20
plus VAT &
Shipping**

Product Images



Description

For vials ND8 with rolled rim.

Transparent PE caps 9 x 5.9 mm with 4 mm hole

Order No. 7-0617: Septum natural rubber red-orange/TEF transparent 1.0 mm thick,

Additional Information

No.	7-0617
Manufacturer (Brand)	neochrom
EAN	4058072028893
Transport temperature	Room temperature
Color	transparent (clear)
Type closure/ lid	Push-on/plug-in closure
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart suitable for autosampler cap with septum
DM Hole	4 mm
Hardness	60 shore A
Material gasket / septum	Natural rubber Polytetrafluoroethylene (PTFE)
Closure / lid properties	with mounted septum with hole
DM Closure	8 mm
Thickness septum	1 mm
Color closure/ lid	transparent (clear)
Material closure/ lid	Low-density polyethylene (LDPE)
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

