



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® PE Push-On Cap blue, ND11 with diluted puncture, 100 pcs/pack

€7.27
plus VAT &
Shipping

Product Images



Description

For vials ND11 with rolled rim.

PE push-on cap with thinned piercing PE, without septum.

Additional Information

No.	7-0652
Manufacturer (Brand)	neochrom
EAN	4058072029173
Transport temperature	Room temperature
Color	blue
Type closure/ lid	Push-on/plug-in closure
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart suitable for autosampler
DM Hole	11 mm
Closure / lid properties	pierceable closed
Color closure/ lid	blue
DM Closure	11 mm
Material closure/ lid	Low-density polyethylene (LDPE)
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

