



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Cryomaster® lid inserts, blue, 500 pcs/pck

€20.90
plus VAT &
Shipping

Product Images



Description

Colorful Cryomaster® lid inserts made of high-quality polypropylene for Cryomaster® cryotubes 4-6100 ff. The lid inserts can be pressed into the lids from above for color coding of the tubes. Thus, a wide range of different color combinations can be made.

36 color combinations by 6 lid colors and 6 insert colors plus the 1D barcode. So you succeed in a fast and clear sample assignment!

PU: 500 pcs.

Additional Information

No.	4-6113
Manufacturer (Brand)	Cryomaster
EAN	4058072154240
Transport temperature	Room temperature
Color	blue
Material	Polypropylene (PP)
suitable for	Closures for cryotubes

