



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Cryomaster® lid inserts, yellow, 500 pcs/pck

**€20.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Colorful Cryomaster® lid inserts made of high-quality polypropylene for Cryomaster® cryotubes 4-6100 ff. The lid inserts can be pressed into the lids from above for color coding of the tubes. Thus, a wide range of different color combinations can be made.

36 color combinations by 6 lid colors and 6 insert colors plus the 1D barcode. So you succeed in a fast and clear sample assignment!

**PU: 500 pcs.**

## Additional Information

---

No.	4-6116
Manufacturer (Brand)	Cryomaster
EAN	4058072154257
Transport temperature	Room temperature
Color	yellow
Material	Polypropylene (PP)
suitable for	Closures for cryotubes

