




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Brand lid inserts for color coding, red; PU: 500 pieces

**€18.75
plus VAT &
Shipping**

Product Images



Description

Lid inserts, red for reaction vessels with attached and loose screw cap. Made of PP. Application range: -196 °C to +121 °C.

Additional Information

No.	BN-2165
Manufacturer (Brand)	Brand
old neoLab article no.	1-6833
EAN	4058072213893
Transport temperature	Room temperature
Color	red
Material	Polypropylene (PP)
suitable for	Closures for reaction vessels

