



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Brand lid inserts for color coding, yellow; PU:
500 pieces**

**€18.75
plus VAT &
Shipping**

Product Images



Description

Lid inserts, yellow for reaction vessels with attached and loose screw cap. Made of PP. Application range: -196 °C to +121 °C.

Additional Information

No.	BN-2167
Manufacturer (Brand)	Brand
old neoLab article no.	1-6835
EAN	4058072213916
Transport temperature	Room temperature
Color	yellow
Material	Polypropylene (PP)
suitable for	Closures for reaction vessels

