



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neochrom® aluminum flare cap (silver), 11 mm , red PTFE/white silicone slotted

€49.70
plus VAT &
Shipping

Product Images



Description

The neochrom aluminum flanged closures ND11 with septum are used to close rolled rim vials ND11.

Additional Information

No.	EC-1024
Manufacturer (Brand)	neochrom
EAN	4058072079130
Transport temperature	Room temperature
Closures and lids features	with slotted septum for blunt. thick needles closure with septum
Material closure/ lid	Aluminum

