



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® aluminum crimp caps ND11 with  
hole, septum silicone white/PTFE blue, 100  
pcs/pack**

**€58.20  
plus VAT &  
Shipping**

## Product Images



## Description

---

### **Aluminum crimp caps ND11, colorless lacquered, with hole 5,5 mm.**

#### **1. with septum natural rubber/TEF**

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.0 mm

#### **2. with septum Red Rubber/PTFE beige**

Temperature resistant from -40°C to +110°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Lower fragmentation than natural rubber/TEF

Thickness 1.0 mm

#### **3. with septum silicone/PTFE**

Temperature resistant from -60°C to +200°C

For single injections

Thickness 1.3 mm

#### **4. with septum PTFE virginal**

High temperature resistant

For single use injections

Thickness 0,25 mm

## Additional Information

No.	EC-1316
Manufacturer (Brand)	neochrom
EAN	4058072086015
Transport temperature	Room temperature
Color	silver
Closures and lids features	with mounted septum suitable for autosampler cap with septum cap autoclavable
DM Hole	5.5 mm
Quantity components	three
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
DM Closure	11 mm
Thickness septum	1.5 mm
Material closure/ lid	Aluminum
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
Type flare cap	Flared cap. standard
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

