



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® flare cap gold ND20 magnetic, hole
8 mm, butyl/PTFE grey, 100 pcs./pack**

**€33.20
plus VAT &
Shipping**

Product Images



Description

Magnetic caps ND20, mounted with septa.

Septum Butyl/PTFE gray

Temperature resistant from -40°C to 120°C

Fully coated with PTFE

Thickness: 3.0 mm

Additional Information

No.	7-0817
Manufacturer (Brand)	neochrom
EAN	4058072030599
Transport temperature	Room temperature
Color	Gold
Type closure/ lid	Flare cap
Closures and lids features	with mounted septum magnetic suitable for autosampler cap with septum cap autoclavable
Quantity components	two
DM Hole	8 mm
DM Closure	20 mm
Thickness septum	3 mm
Hardness	55 shore A
Material gasket / septum	Polytetrafluoroethylene (PTFE)
Material closure/ lid	Steel
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
Type flare cap	Flare cap. magnetic
Closure / lid properties	with hole with mounted septum

