



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® Alu-flange cap ND20, mill finish,  
hole 10 mm, Pharma-Fix Butyl/ PTFE, 100  
pcs./pack**

**€41.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

Flare closures, assembled with septa.

Septum Butyl/PTFE gray

Temperature resistant from -40°C to 120°C

Fully coated with PTFE

Thickness: 3.0 mm

Pharma-Fix septa

Special shaped discs coated with PTFE only in the contact area of the specimen.

On the glass edges, the elastic butyl or silicone achieves very good sealing.

Temperature resistance see above

Thickness: 3.0 mm

## Additional Information

No.	7-0820
Manufacturer (Brand)	neochrom
EAN	4058072030629
Transport temperature	Room temperature
Color	silver
Type closure/ lid	Flare cap
Closures and lids features	with mounted septum suitable for autosampler cap with septum cap autoclavable septum especially tight sealing (Pharma-Fix)
DM Hole	10 mm
Quantity components	two
Hardness	55 shore A
Material gasket / septum	Polytetrafluoroethylene (PTFE)
Closure / lid properties	with hole with mounted septum
DM Closure	20 mm
Thickness septum	3 mm
Material closure/ lid	Aluminum
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
Type flare cap	Flared cap. standard
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Agilent G1888A Agilent HS 7694 Carlo Erba HS850 Dani HS 39.50 Dani HS 86.50 Fisons HS850 Perkin Elmer HS6 Thermo Scientific HS850 Thermo Scientific HS2000 Thermo Scientific TriPlus HS Teledyne Tekmar HT3 Varian CP-9020/9025 Varian CP-9060 Varian Genesis

