



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Magnet. Crimp cap ND20, gold with
hole, silicone blue/PTFE transp., 100 pcs./pack**

**€59.20
plus VAT &
Shipping**

Product Images



Description

Magnetic caps ND20, mounted with septa.

Septum Silicone/PTFE

Temperature resistant from -60°C to 200°C

For easy penetration

For critical analysis

Thickness: 3.0 mm

Additional Information

No.	7-0827
Manufacturer (Brand)	neochrom
EAN	4058072030681
Transport temperature	Room temperature
Color	Gold
Type closure/ lid	Flare cap
Closures and lids features	with mounted septum with high purity septum for critical analyses magnetic suitable for autosampler cap with septum cap autoclavable
Quantity components	two
DM Hole	5 mm
DM Closure	20 mm
Thickness septum	3 mm
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Material closure/ lid	Steel
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Type flare cap	Flare cap. magnetic
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Carlo Erba HS250 Carlo Erba HS500 Carlo Erba HS800 CTC HS 500 Fisons HS250 Fisons HS500 Fisons HS800 Thermo Scientific HS250 Thermo Scientific HS500 Thermo Scientific HS800

