



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® Headspace cap ND 20 colorless,  
overpressure protection, silicone white/PTFE  
beige, 100 pcs./pack**

**€58.60  
plus VAT &  
Shipping**

## Product Images



## Description

---

Flare closures, assembled with septa.  
With overpressure protection/set break point.

Septum silicone/PTFE

Temperature resistant from -60°C to 200°C

For easy penetration

For critical analyses

Thickness: 3.0 mm

Silicone white/PTFE beige = AGILENT quality, thickness 3.2 mm

## Additional Information

No.	7-0832
Manufacturer (Brand)	neochrom
EAN	4058072030711
Transport temperature	Room temperature
Color	silver
Type closure/ lid	Flare cap
Closures and lids features	with overpressure protection with mounted septum suitable for autosampler suitable for headspace applications cap with septum cap autoclavable
Quantity components	two
DM Hole	8 mm
DM Closure	20 mm
Thickness septum	3.2 mm
Hardness	45 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Material closure/ lid	Aluminum
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Type flare cap	Headspace cap
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Agilent G1888A Agilent HS 7694 Dani HS 39.50 Dani HS 86.50 Perkin Elmer HS6

