



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Bohlender disposable ground joint sleeves (PTFE) NS 45/40, 5 pcs./pack

€41.55 As low
as
€31.85
plus VAT &
Shipping



Description

Sealing without grease!

Excellent chemical resistance at service temperature from -200°C to +260°C.

Wall thickness: approx. 0.05 mm

Vacuum: suitable up to 0.1 bar

Additional Information

No.	ER-1658
Manufacturer (Brand)	Bohlander
old neoLab article no.	1-5149
EAN	4058072245450
Transport temperature	Room temperature
Material	Polytetrafluoroethylene (PTFE)
Thickness	0.05 mm
TBST MAX	260 °C
for NS	45/40

