



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® screw cap ND9 short thread with
hole, natural rubber red-or./TEF transp. 100
pcs./pack**

**€17.10
plus VAT &
Shipping**

Product Images



Description

Screw cap ND9 made of PP, with short thread 9-425 and hole 6 mm.

With septum natural rubber/TEF

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.0 mm

Additional Information

No.	7-0689
Manufacturer (Brand)	neochrom
EAN	4058072029524
Transport temperature	Room temperature
Color	transparent (milky)
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart suitable for autosampler cap with septum cap autoclavable
DM Hole	6 mm
Thread	Short thread
Hardness	60 shore A
Material gasket / septum	Natural rubber Polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	with hole with mounted septum
DM Closure	9 mm
Thickness septum	1 mm
Color closure/ lid	transparent (milky)
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

