



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® short thread screw cap ND9, PP red  
with hole, rubber red-or./TEF transp. 100  
pcs./pack**

**€17.70  
plus VAT &  
Shipping**

## Product Images



## Description

---

Red screw cap ND9 made of PP, with short thread 9-425 and hole 6 mm.

With septum natural rubber/TEF

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.0 mm

## Additional Information

No.	7-0691
Manufacturer (Brand)	neochrom
EAN	4058072029548
Transport temperature	Room temperature
Color	red
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart with slotted septum for blunt, thick needles suitable for autosampler suitable for multiple injections cap with septum cap autoclavable
DM Hole	6 mm
Thread	Short thread
Quantity components	two
Hardness	60 shore A
Material gasket / septum	Natural rubber Polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	with hole with mounted septum
DM Closure	9 mm
Thickness septum	1 mm
Color closure/ lid	red
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

