




neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Screw cap with seal, transparent, 500  
pcs./pack**

**€80.00  
plus VAT &  
Shipping**

## Product Images



## Description

---

**Transparent screw caps with sealing ring** that fit any reaction vessel.

The lid is slightly ribbed for ease of work. The O-ring is made of EPR (ethylene propylene).

Packaging: 500 pieces in a bag.

## Additional Information

---

No.	1-6193
Manufacturer (Brand)	neoLab
EAN	4058072042844
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	with gasket
Color closure/ lid	transparent (clear)
Material gasket / septum	Ethylene propylene rubber (EPR)
Closure / lid properties	closed

