



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Screw cap with seal, orange, 500 pcs./pack

**€80.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**Orange colored screw cap with sealing ring** that fits on any reaction vessel.

The lid is slightly ribbed for ease of work. The O-ring is made of EPR (ethylene propylene).

Packaging: 500 pieces in a bag.

## Additional Information

---

No.	1-6194
Manufacturer (Brand)	neoLab
EAN	4058072042851
autoclavable at 121°C	Yes
Color	orange
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	with gasket
Material gasket / septum	Ethylene propylene rubber (EPR)
Closure / lid properties	closed
Color closure/ lid	orange

