




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Screw cap with seal, yellow, 500 pcs./pack

€80.00
plus VAT &
Shipping

Product Images



Description

Yellow screw caps with sealing ring that fit any reaction vessel.

The lid is slightly ribbed for ease of work. The O-ring is made of EPR (ethylene propylene).

Packaging: 500 pieces in a bag.

Additional Information

No.	1-6197
Manufacturer (Brand)	neoLab
EAN	4058072042882
autoclavable at 121°C	Yes
Color	yellow
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	with gasket
Material gasket / septum	Ethylene propylene rubber (EPR)
Closure / lid properties	closed
Color closure/ lid	yellow

