



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine screw cap GL 45 with 2 integrated hose connections

€16.90
plus VAT &
Shipping

Product Images



Description

Screw cap GL 45 with silicone seal and 2 integrated hose connections.

(1 x air/1 x liquid)

This allows any laboratory bottle with the appropriate thread size (possibly Ex-protected) to be converted into a suction bottle.

Additional Information

No.	2-9346
Manufacturer (Brand)	neoLabLine
EAN	4058072086770
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	blue
Material	Polypropylene (PP)
sterile	No
Thread	GL 45
Number of valves / cocks	2
for DM hose inside 1	6-10 mm
DM Closure	45 mm
Color closure/ lid	blue
Type connection collection bottle	Hose olive

