



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Closing cap GL 45, with pouring ring, PP, yellow, PU 10 pieces

**€16.20**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Screw caps GL 45 made of PP with lip seal and pouring ring.  
In yellow.  
Steam sterilizable up to 121°C.

## Additional Information

---

No.	2-3065
Manufacturer (Brand)	neoLab
EAN	4058072000516
Transport temperature	Room temperature
Color	yellow
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	with seal cap autoclavable
Thread	GL 45
Closure / lid properties	with seal closed
DM Closure	45 mm
Color closure/ lid	yellow
for bottle shape	round eng

