



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Closing cap GL 45, with pouring ring, PP, purple, pack of 10

€16.20
plus VAT &
Shipping

Product Images



Description

Screw caps GL 45 made of PP with lip seal and pouring ring.
In purple.
Steam sterilizable up to 121°C.

Additional Information

No.	2-3067
Manufacturer (Brand)	neoLab
EAN	4058072000530
Transport temperature	Room temperature
Color	purple
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	with seal cap autoclavable
Thread	GL 45
Closure / lid properties	with seal closed
DM Closure	45 mm
Color closure/ lid	purple
for bottle shape	round eng

