



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Screw Cap Poly Insert DIN 55

€0.32 As low as

€0.29

**plus VAT &
Shipping**

Product Images



Description

Black screw thread caps with plastic insert,

for wide neck glass bottles .

Accessories for: Wide-neck bottles (Order No. 2-1471 ff, 2-1480 ff).

Additional Information

No.	2-1497
Manufacturer (Brand)	neoLab
EAN	4058072057039
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Height	17.5 mm
Type closure/ lid	Screw cap
DM outside	59 mm
Closures and lids features	Closure with septum
Thread	DIN 55
Closure / lid properties	closed
DM Closure	55 mm
Color closure/ lid	black
for bottle shape	round wide

