



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® screw cap GL 32 made of PP blue, with lip seal, 10 pcs./pack

€17.30
plus VAT &
Shipping

Product Images



Description

As accessories for ^{DURAN®} laboratory bottles we offer:

Screw caps with lip seal.

Made of blue PP, autoclavable and temperature resistant up to +140°C.

Additional Information

| | |
|----------------------------|----------------------------|
| No. | SD-2748 |
| Manufacturer (Brand) | DURAN |
| old neolab article no. | 2-1393 |
| EAN | 4058072363543 |
| Transport temperature | Room temperature |
| Color | blue |
| Material | Polypropylene (PP) |
| Height | 24 mm |
| Type closure/ lid | Screw cap |
| DM outside | 40 mm |
| Closures and lids features | with seal cap autoclavable |
| DM Closure | 32 mm |
| Color closure/ lid | blue |
| Thread | GL 32 |
| Closure / lid properties | with seal closed |
| for bottle shape | round eng |

