



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® screw cap GL 45 made of PP blue, with lip seal, 10 pcs./pack

€20.60
plus VAT &
Shipping

Product Images



Description

As accessories for ^{DURAN®} laboratory bottles we offer:

Screw caps with lip seal.

Made of blue PP, autoclavable and temperature resistant up to +140°C.

Additional Information

No.	SD-2749
Manufacturer (Brand)	DURAN
old neoLab article no.	2-1394
EAN	4058072363550
Transport temperature	Room temperature
Color	blue
Material	Polypropylene (PP)
Height	25 mm
Type closure/ lid	Screw cap
DM outside	54 mm
Closures and lids features	with seal cap autoclavable
Thread	GL 45
Closure / lid properties	with seal closed
Color closure/ lid	blue
DM Closure	45 mm
for bottle shape	round eng

