



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Screw caps made of PPN blue, GL 45 10 pcs./pack blue

**€13.80**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Screw caps suitable for<sup>Duran®</sup> pressure plus laboratory bottle with DIN thread GL 45

**Order No. E-1433**

**Material:** PPN/blue  
closed cap, max. + 140°C

## Additional Information

---

No.	E-1433
Manufacturer (Brand)	neoLab
EAN	4058072064730
Transport temperature	Room temperature
Color	blue
Material	Polypropylene (PP)
Height	25 mm
Type closure/ lid	Screw cap
DM outside	54 mm
Closures and lids features	with seal cap autoclavable
Thread	GL 45
Closure / lid properties	closed
Color closure/ lid	blue
DM Closure	45 mm
for bottle shape	round eng

