



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® screw cap GL 45 with hole red up to +180 °C

€10.66 As low as
€8.20
plus VAT &
Shipping

Product Images



Description

Screw caps suitable for^{Duran®} pressure plus laboratory bottle with DIN thread GL 45

Order No. 2-3051

Material: PBT/red

with bore Ø 34 mm, max.+ 180°C

Additional Information

No.	SD-2729
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3051
EAN	4058072363369
Transport temperature	Room temperature
Color	red
Material	Polybutylene terephthalate (PBT)
Height	25 mm
Type closure/ lid	Screw cap
DM outside	54 mm
Closures and lids features	with seal cap autoclavable
Thread	GL 45
Material gasket / septum	Perfluoroethylene propylene (FEP) Silicone
Closure / lid properties	with hole
Color closure/ lid	red
DM Closure	45 mm
for bottle shape	round eng

