



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® screw cap magnetic ND18 silver,
8mm hole, silicone blue/PTFE white (pack á 100
pieces)**

**€70.20
plus VAT &
Shipping**

Product Images



Description

Magnetic screw caps ND18, silver, mounted with septa.

With hole 8 mm Ø

With septum silicone blue/PTFE white

Temperature resistant from -60°C to 200°C

Additional Information

| | |
|--|---|
| No. | 7-0845 |
| Manufacturer (Brand) | neochrom |
| EAN | 4058072037130 |
| Transport temperature | Room temperature |
| Color | silver |
| Material | Steel |
| Type closure/ lid | Screw cap |
| Closures and lids features | Reference to non-verified autosampler compatibility in shopping cart with high purity septum for critical analysis magnetic suitable for autosampler suitable for headspace applications cap with septum cap autoclavable |
| DM Hole | 8 mm |
| Thread | Fine thread |
| Quantity components | two |
| Hardness | 45 shore A |
| Material gasket / septum | Silicone Polytetrafluoroethylene (PTFE) |
| Texture | corrugated |
| Closure / lid properties | with hole with mounted septum |
| DM Closure | 18 mm |
| Thickness septum | 1.3 mm |
| Color closure/ lid | silver |
| TBST MAX Septum | 200 °C |
| TBST MIN Septum | -60 °C |
| for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS) | CTC (LEAP) Combi PAL/ GC PAL. 32 pos. tray CTC (LEAP) Combi PAL/ GC PAL. 98 Pos. Tray CTC (LEAP) Combi PAL/ GC PAL. 200 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 32 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 98 Pos. Tray CTC Combi-xt Headspace Option CTC C |

