



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® screw cap magnetic ND18 silver,
8mm hole, silicone white/PTFE blue (pack á 100
pieces)**

**€86.90
plus VAT &
Shipping**

Product Images



Description

Magnetic screw caps ND18, silver, mounted with septa.
With hole 8 mm Ø

With septum silicone white/PTFE blue
Temperature resistant from -60°C to 200°C

Additional Information

No.	7-0846
Manufacturer (Brand)	neochrom
EAN	4058072037147
Transport temperature	Room temperature
Color	silver
Material	Steel
Type closure/ lid	Screw cap
Closures and lids features	Reference to non-verified autosampler compatibility in shopping cart with high purity septum for critical analysis magnetic suitable for autosampler suitable for headspace applications cap with septum cap autoclavable
DM Hole	8 mm
Thread	Fine thread
Quantity components	two
Hardness	55 shore A
Material gasket / septum	Silicone Polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	with hole with mounted septum
DM Closure	18 mm
Thickness septum	1.5 mm
Color closure/ lid	silver
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Agilent G1888A Agilent CTC Combi PAL CTC (LEAP) Combi PAL/ GC PAL. 32 pos. tray CTC (LEAP) Combi PAL/ GC PAL. 98 Pos. Tray CTC (LEAP) Combi PAL/ GC PAL. 200 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 32 Pos. Tray CTC (LEAP) Combi PAL SMPE Mode. 98 Pos. Tra

