



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® screw cap magnetic ND18 silver,
8mm hole, butyl red/PTFE grey (pack á 100
pieces)**

**€40.20
plus VAT &
Shipping**

Product Images



Description

Magnetic screw caps ND18, silver, mounted with septa.
With hole 8 mm Ø

With septum Butyl red/PTFE grey
Temperature resistant from -40°C to 120°C

Additional Information

No.	7-0849
Manufacturer (Brand)	neochrom
EAN	4058072037178
Transport temperature	Room temperature
Color	silver
Material	Steel
Type closure/ lid	Screw cap
Closures and lids features	Note on non-verified autosampler compatibility in shopping cart magnetic suitable for autosampler suitable for headspace applications cap with septum cap autoclavable
DM Hole	8 mm
Thread	Fine thread
Quantity components	two
Hardness	55 shore A
Material gasket / septum	Butyl rubber (IIR) polytetrafluoroethylene (PTFE)
Texture	corrugated
Closure / lid properties	with hole with mounted septum
DM Closure	18 mm
Thickness septum	1.6 mm
Color closure/ lid	silver
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Perkin Elmer Turbomatrix HS16/ HS40/ HS40 XL/ HS40 Trap/ HS110/ HS110 Trap

