



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® Magn. Screw caps ND18, silver, hole
8 mm, silicone white/aluminum foil silver, 100
pcs./pack**

**€113.10
plus VAT &
Shipping**

Product Images



Description

Magnetic screw caps ND18, silver, mounted with septa.

With hole 8 mm Ø

1. with septum silicone/PTFE

Temperature resistant from -60°C to 200°C

2. with septum Butyl/PTFE

Temperature resistant from -40°C to 120°C

Additional Information

No.	EC-1326
Manufacturer (Brand)	neochrom
EAN	4058072086114
Transport temperature	Room temperature
Color	silver
Material	Steel
Closures and lids features	magnetic suitable for autosampler suitable for headspace applications cap with septum cap autoclavable
Quantity components	two
Texture	corrugated
DM Hole	8 mm
DM Closure	18 mm
Thickness septum	1.3 mm
Color closure / lid	silver
Thread	Fine thread
Hardness	50 shore A
Material gasket / septum	Silicone Aluminum foil
TBST MAX Septum	200 °C
TBST MIN Septum	-60 °C
Closure / lid properties	closed with mounted septum

