



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neochrom® screw caps PP, ND10 w. Hole, natural  
rubber red-orange / TEF transparent Thickness  
1,3 mm, 60° shore A 100 pcs./pack**

**€19.30  
plus VAT &  
Shipping**

## Product Images



## Description

---

Screw caps ND10 made of PP, with thread 10-425 and hole 7 mm.

With septum natural rubber red-orange/TEF transparent

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to the good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

## Additional Information

No.	7-0728
Manufacturer (Brand)	neochrom
EAN	4058072029821
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Type closure/ lid	Screw cap
Closures and lids features	Reference to non-verified autosampler compatibility in shopping cart suitable for autosampler suitable for multiple injections cap with septum cap autoclavable
Quantity components	two
Texture	corrugated
DM Hole	7 mm
DM Closure	10 mm
Thickness septum	1.3 mm
Color closure/ lid	black
Thread	10-425 (US)
Hardness	60 shore A
Material gasket / septum	Natural rubber Polytetrafluoroethylene (PTFE)
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
Closure / lid properties	with hole with mounted septum
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

