



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neochrom® screw caps PP, ND10 closed, natural
rubber red-orangeSilic. white/TEF transparent,
100 pcs./pack**

**€21.10
plus VAT &
Shipping**

Product Images



Description

Screw caps ND10 made of PP, with thread 10-425 and hole 7 mm.

1. with septum natural rubber/TEF

Temperature resistant from -40°C to +120°C

Ideal for multiple injections due to good resealing properties

Standard closure for GC and HPLC

Thickness 1.3 mm

2. with septum silicone/PTFE

Temperature resistant from -60°C to +200°C

For single injections

Additional Information

No.	EC-1309
Manufacturer (Brand)	neochrom
EAN	4058072085995
Transport temperature	Room temperature
Color	Black
Material	Polypropylene (PP)
Closures and lids features	suitable for autosampler suitable for multiple injections cap with septum cap autoclavable
Thread	10-425 (US)
Quantity components	two
Hardness	60 shore A
Material gasket / septum	Natural rubber
Texture	corrugated
Closure / lid properties	closed with mounted septum
Color closure/ lid	black
DM Closure	10 mm
Thickness septum	1.3 mm
TBST MAX Septum	120 °C
TBST MIN Septum	-40 °C
for the autosampler model (If your AS does not appear here, please ask the AS manufacturer directly for the compatibility of the cap with your AS)	Please clarify the compatibility of this cap with your autosampler directly with the manufacturer of your autosampler

